

ANALYSIS OF VINEGAR
Experiment #5


PURPOSE:

- To determine the percent by weight of acetic acid in vinegar.
- To perform an Acid-Base titration.

PRINCIPLES:

Vinegar is an aqueous solution of acetic acid ($\text{HC}_2\text{H}_3\text{O}_2$) commonly used in flavoring and preserving food. This is because the acidity adds “tang” and inhibits bacterial growth.

In the investigation that follows the strength of one brand of vinegar will be determined by titration. The acetic acid is neutralized with a standardized sodium hydroxide solution as shown in the following equation:



Any vinegar sample may be used, but colorless vinegar is preferred because it gives less interference with the observation of the indicator endpoint color change.

As the concentration of the vinegar solution is much higher than the concentration of your standardized sodium hydroxide solution, the original vinegar solution is diluted ten times prior to the titration. This dilution factor must be taken into account when calculating the concentration of the original vinegar solution.

PROCEDURE:

1. Clean and dry the following material:
 - a. 500 mL plastic bottle ❶
 - b. 100 mL beaker ❷
2. Clean and rinse the following glassware with distilled water
 - a. 25.00 ml volumetric pipet ❸
 - b. 250.0 mL volumetric flask ❹ **(CAUTION: These are expensive!!)**
 - c. 250 mL Erlenmeyer flask ❺
 - d. 50 mL Buret ❻

3. Measure exactly 25.00 mL ③ of vinegar into a clean 250.0 mL volumetric flask ④.

Notes:

- Be sure to use a dry beaker to transfer the vinegar.
 - Do not put any used vinegar back in the vinegar supply.
 - Be sure to rinse your pipette with vinegar prior to using it.
4. Dilute the vinegar with deionized water to the mark on the volumetric flask ④.
 5. Stopper the flask ④ and mix the solution well. (Invert the solution slowly for at least 10 times to completely mix the contents).
 6. Transfer the dilute vinegar solution to a clean and dry 500 mL plastic bottle ① and label it with both the contents and your name.
 7. Immediately wash your volumetric flask ④ with plenty of tap water and several portions of deionized water. Let the flask dry at room temperature.
 8. Pour about 50 mL of the dilute vinegar solution in your 100 mL beaker ②.
 9. Rinse your 25.00 mL volumetric pipet ③ several times with portions of diluted vinegar from your beaker.



Be careful not accidentally add any water to the diluted vinegar solution in the beaker.



10. Carefully pipet 25.00 mL ③ of diluted vinegar solution into the 250 mL Erlenmeyer flask ⑤.
11. Add about 50 mL of deionized water to the Erlenmeyer flask ⑤.
12. Add 2 drops of phenolphthalein indicator solution and swirl the flask to thoroughly mix the solution.
13. Rinse your 50 mL buret ⑥ several times with a few milliliters of your standardized sodium hydroxide solution.
14. Fill the buret ⑥ with your standardized sodium hydroxide solution.
 - a. Make sure that the tip does not have any air bubbles.
 - b. Record the volume on the buret to the nearest 0.01 mL.
15. Titrate the acid sample to a faint pink end point.
16. Record the final volume of the buret to the nearest 0.01 mL.

17. Repeat the titration procedure described above for at least two more trials. The number of trials run depends on:
- How much standardized sodium hydroxide solution you have available
 - The precision of your data.

CALCULATIONS:

$$\text{Moles of base} = (\text{Molarity of base}) \times (\text{Liters of base})$$

$$\text{Moles of acid (diluted vinegar)} = \text{Moles of base}$$

$$\text{Molarity of Diluted Vinegar} = \frac{\text{Moles of acid (diluted vinegar)}}{\text{Liters of acid (diluted vinegar)}}$$

The **MOLARITY OF THE ORIGINAL VINEGAR SOLUTION** can be found by keeping in mind that the vinegar has been diluted ten times (from 25.00 mL to 250.0 mL) to obtain the **DILUTED VINEGAR**, whose molarity has been determined by titration with standardized NaOH.

$$\frac{\text{Grams of } HC_2H_3O_2}{\text{Grams of vinegar}} = \frac{\text{mol } HC_2H_3O_2}{L \text{ of vinegar}} \times \frac{g \text{ } HC_2H_3O_2}{\text{mol } HC_2H_3O_2} \times \frac{L \text{ vinegar}}{mL \text{ vinegar}} \times \frac{mL \text{ vinegar}}{g \text{ vinegar}}$$

$$\text{Density of vinegar} = 1.0052 \text{ g / mL}$$

$$\frac{\% \text{ Weight } HC_2H_3O_2}{\text{Weight vinegar}} = \frac{\text{Grams of } HC_2H_3O_2}{\text{Grams of vinegar}} \times 100$$

EXPERIMENT #5
REPORT FORM

Molarity of Standardized NaOH solution: _____ M

	TITRATIONS		
	First	Second	Third
Final buret reading (mL)	_____	_____	_____
Initial buret reading (mL)	_____	_____	_____
Volume of titrant (mL)	_____	_____	_____
Moles NaOH used (mol)	_____	_____	_____
<u>Show one sample calculations below:</u>			

Moles of HC ₂ H ₃ O ₂ used	_____	_____	_____
Volume of acid used (L) (diluted vinegar)	_____	_____	_____
Molarity of acid (M) (diluted vinegar)	_____	_____	_____
<u>Show one sample calculations below:</u>			

Average molarity of acid (M) (diluted vinegar)	_____
Deviation d _i	_____
Mean deviation	_____
Relative Mean Deviation, RMD	_____
Molarity of original vinegar (M) (before dilution)	_____
$\frac{\text{Weight of HC}_2\text{H}_3\text{O}_2}{\text{Weight of vinegar}}, \%$	_____